

mediven[®] struva[®] 23

Pressure of 23 mmHg in the ankle area



mediven struva 23

Wound dressing aids have a comprehensive range of uses in the outpatient and partial inpatient treatment of thrombosis, in vein surgery and as a general replacement for compression bandages for the lower extremities.

With a pressure of 23 mmHg in the ankle area, which decreases continuously towards the thigh, mediven struva 23 is particularly suitable for postoperative short term compression directly after vein surgery. Venous return is promoted and the development of oedema and effusions prevented.

In vein surgery, mediven struva 23 is used as a supplement to mediven struva 35

- During and after the operation, the unoperated leg is fitted with the mediven struva 23 to prevent thrombosis.
- After one to four days, compression treatment of the operated leg with mediven struva 35 is discontinued and this can then also be fitted with a mediven struva 23.

Material: 76 % polyamide, 24 % elastane, latex-free